




## CORDS

Since 1979, Kirwan Surgical Products has built its reputation on manufacturing superior quality electrosurgical accessories.

Kirwan Surgical cords, renowned for their flexibility, suppleness and durability, are unsurpassed in the industry.

| ITEM #  | LENGTH     | DESCRIPTION  |  |
|---------|------------|--|--|
| 10-6010 | 12' (3.6m) | Reusable silicone bipolar cord with Lemo Plug (for Storz-type Daisy)                         |  |
| 10-6011 | 12' (3.6m) | Reusable silicone bipolar cord with Lemo Plug (for Storz-type Premier)                       |  |
| 10-6012 | 12' (3.6m) | Reusable silicone bipolar cord with 2 Lemo plugs   |  |
| 10-6020 | 12' (3.6m) | Reusable silicone bipolar cord (with banana pins & bladed connector)                         |  |
| 10-6025 | 12' (3.6m) | Reusable silicone bipolar cord with booted forceps connector (for Wolf/Martin-type)          |  |
| 10-6050 | 12' (3.6m) | Reusable silicone bipolar cord with Erbe-type/Storz Millennium-type booted forceps connector |  |
| 10-6055 | 12' (3.6m) | Reusable silicone bipolar cord with Redel connector  |  |
| 10-6060 | 12' (3.6m) | Reusable silicone bipolar cord with Fischer plug   |  |



## REUSABLE BIPOLAR CORDS

| ITEM #  | LENGTH      | DESCRIPTION   |
|---------|-------------|---|
| 10-6075 | 12' (3.6m)  | Reusable silicone bipolar cord (for O.M.S. -type Phaco machine) |
| 10-6076 | 6' (1.8m)   | Reusable silicone bipolar cord (for O.M.S. -type Phaco machine) |
| 10-7001 | 12' (3.6m)  | Reusable silicone bipolar handswitch cord                       |
| 10-7050 | 6'4" (1.9m) | Reusable blue silicone coaxial cord (for D.O.R.C. -type Phaco)  |







## DISPOSABLE FIXED PIN BIPOLAR CORDS

| ITEM #   | LENGTH     | DESCRIPTION  | FEATURES |
|----------|------------|--|----------|
| 10-4000J | 12' (3.6m) | Disposable bipolar cord with fixed pins, 38mm, sterile green   |          |
| 10-4001J | 12' (3.6m) | Disposable bipolar cord with fixed pins, 38mm, sterile clear   |          |
| 10-4000R | 12' (3.6m) | Disposable bipolar cord with fixed pins, 22mm, sterile green   |          |
| 10-4001R | 12' (3.6m) | Disposable bipolar cord with fixed pins, 22mm, sterile clear   |          |
| 10-4000V | 12' (3.6m) | Disposable bipolar cord with fixed pins, 28.6mm, sterile green |          |
| 10-4001V | 12' (3.6m) | Disposable bipolar cord with fixed pins, 28.6mm, sterile clear |          |



## DISPOSABLE FIXED PIN BIPOLAR CORDS

| ITEM #   | LENGTH | DESCRIPTION   | FEATURES  |   |
|----------|--------|---|---|---|
| 10-4601J |        | Disposable bipolar irrigation tube and cord assembly with fixed pins, 38mm, sterile             |  |  |
| 10-4602J |        | Disposable bipolar irrigation tube and cord assembly, M1000 type with fixed pins, 38mm, sterile |  |   |
| 10-4601V |        | Disposable bipolar irrigation tube and cord assembly with fixed pins, 28.6mm, sterile           |  |   |






## REUSABLE FIXED PIN BIPOLAR CORDS

| ITEM #   | LENGTH     | DESCRIPTION  |   |
|----------|------------|--|---|
| 10-5000R | 12' (3.6m) | Reusable TPR bipolar cord with fixed pins, 22mm                    |  |
| 10-6000R | 12' (3.6m) | Reusable silicone bipolar cord with fixed pins, 22mm               |   |
| 10-5000V | 12' (3.6m) | Reusable TPR bipolar cord with fixed pins, 28.6mm                  |  |
| 10-6000V | 12' (3.6m) | Reusable silicone bipolar cord with fixed pins, 28.6mm             |   |
| 10-6000J | 12' (3.6m) | Reusable silicone bipolar cord with fixed pins, 38mm               |  |
| 10-6020J | 12' (3.6m) | Reusable silicone bipolar cord with fixed pins, 38mm, EU Connector |  |





# CORD ADAPTERS

| ITEM #  | DESCRIPTION  |  |
|---------|--|--|
| 10-6500 | Reusable bipolar Erbe-type/Storz-Millennium-type adapter (Erbe-type to banana pin) |   |
| 10-6575 | Reusable O.M.S.-type Phaco machine adapter   |   |
| 30-5110 | Reusable monopolar suction coagulator adapter (banana pin to Bovie)                |   |
| 30-5205 | Reusable monopolar suction coagulator adapter (Resectoscope to Bovie)              |  |
| 30-5210 | Reusable monopolar suction coagulator adapter (Resectoscope to banana pin)         |  |

# DISPOSABLE MONOPOLAR CORDS

| ITEM #   | LENGTH   | DESCRIPTION  | FEATURES  |   |
|----------|----------|--|---|---|
| 30-1201  | 10' (3m) | Disposable monopolar endoscopic cord, sterile green with straight Kirwan plug          |  |  |
| 30-1201L | 10' (3m) | Disposable monopolar endoscopic cord, sterile green with extended straight Kirwan plug |  |   |
| 30-1203  | 10' (3m) | Disposable monopolar endoscopic cord, sterile green with Kirwan plug                   |  |   |
| 30-1204  | 10' (3m) | Disposable monopolar endoscopic cord, sterile green with Kirwan plug and banana pin    |  |   |




## DISPOSABLE MONOPOLAR CORDS

| ITEM #  | LENGTH   | DESCRIPTION   | FEATURES  |   |
|---------|----------|---|---|---|
| 30-1208 | 10' (3m) | Disposable monopolar endoscopic cord, sterile green with Luer connector |  |  |
| 30-1209 | 10' (3m) | Disposable monopolar endoscopic cord, sterile green with 2mm pin        |  |  |

## REUSABLE MONOPOLAR CORDS

| ITEM #  | LENGTH   | DESCRIPTION   |   |
|---------|----------|---|---|
| 30-1000 | 10' (3m) | Reusable silicone monopolar cord with Kirwan Plug                           |  |
| 30-1205 | 10' (3m) | Reusable silicone monopolar endoscopic cord with Kirwan Plug                |  |
| 30-1206 | 10' (3m) | Reusable silicone monopolar endoscopic cord with Kirwan Plug and banana pin |  |
| 30-3003 | 10' (3m) | Reusable silicone monopolar cord with banana pins                           |  |

# INTERNATIONAL BIPOLAR CORDS

| ITEM #   | DESCRIPTION   | FEATURES  |   |
|----------|---|---|---|
| 10-4018J | Disposable 5m bipolar cord, sterile, fixed 38mm pins to European bladed connector   |    |    |
| 10-4021V | Disposable 5m bipolar cord, sterile, fixed 28.6mm pins to European bladed connector |    |    |
| 10-4028J | Disposable 5m bipolar cord, sterile, fixed 38mm pins to booted dual connector       |    |    |
| 10-4028V | Disposable 5m bipolar cord, sterile, fixed 28.6mm pins to booted dual connector     |    |   |
| 10-4038  | Disposable 5m bipolar cord, sterile, Erbe-type to European bladed connector         |    |   |
| 10-4048  | Disposable 5m bipolar cord, sterile, Wolf/Martin-type to European bladed connector  |  |  |
| 10-4058  | Disposable 5m bipolar cord, sterile, Erbe-type to booted dual connector             |  |  |
| 10-4068  | Disposable 5m bipolar cord, sterile, Wolf/Martin-type to booted dual connector      |  |  |