



Access[®] Arthroscopy instruments are designed with more than 50 different tip styles and multiple shaft configurations

- Access to the anterior, middle and posterior areas of the meniscus
- Ergonomic loop handle designed with the surgeon's comfort in mind
- Can be used with either finger loops or palmed for comfort and a balanced feel



Access[®] Knee Arthroscopy Instruments



Linear Punches

- Our linear punches are manufactured with a patented pinned-hinge design for maximum strength and durability.
- The open slot design aids in effective cleaning and sterilization.
- Punches with a meniscal list promote increased access to the lower jaw, under meniscal rim and allow easier cutting of tissue.

Probes

CODE	DESCRIPTION	WORKING DIAMETER
1010501	Probe, Straight, 90° Tip	4.9mm
1010502*	Probe, Curved Left 30°	8.8mm
1010503	Probe, Curved Right 30°	8.8mm



Linear Punches


CODE	DESCRIPTION	BITE WIDTH	WORKING DIAMETER
1010101*	1.7mm, Straight	1.7mm	4.4mm
1010102	1.7mm, Up 15°	1.7mm	4.4mm
1010103	1.7mm, Straight, Meniscal Lift	1.7mm	4.4mm
1010104	1.7mm, Up 15°, Meniscal Lift	1.7mm	4.4mm



CODE	DESCRIPTION	BITE WIDTH	WORKING DIAMETER
1010113	3.6mm, Straight	3.6mm	5.2mm
1010114*	3.6mm, Straight, Meniscal Lift	3.6mm	5.2mm
1010115	3.6mm, Up 15°	3.6mm	5.2mm
1010116	3.6mm, Up 15°, Meniscal Lift	3.6mm	5.2mm



CODE	DESCRIPTION	BITE WIDTH	WORKING DIAMETER
1010120	0.8mm, Straight	0.8mm	4.4mm
1010121*	0.8mm, Up 15°	0.8mm	4.4mm



CODE	DESCRIPTION	BITE WIDTH	WORKING DIAMETER
1010122	2.4mm, Straight, Retro Left	2.4mm	6.6mm
1010123	2.4mm, Straight, Retro Right	2.4mm	6.6mm




Linear Punches

CODE	DESCRIPTION	BITE WIDTH	WORKING DIAMETER
1010112	1.7mm, Up 15°, Blunt Nose	1.7mm	4.4mm



CODE	DESCRIPTION	BITE WIDTH	WORKING DIAMETER
1010133	5.0mm, Straight, Round	5.0mm	4.4mm



CODE	DESCRIPTION	BITE WIDTH	WORKING DIAMETER
1010135	1.7mm, Up 15°, Short Nose	1.7mm	4.4mm



90° Rotary Punches

CODE	DESCRIPTION	BITE WIDTH	WORKING DIAMETER
1010124	90° Tip, Left	3.2mm	7.0mm
1010127	90° Tip, Right	3.2mm	7.0mm




Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use.

Visit symmetrysurgical.com or contact your local Sales Representative for more information.

Rotary Scissors

CODE	DESCRIPTION	BITE WIDTH	WORKING DIAMETER
1010202	20° Tip, Hooked Left	N/A	6.8mm
1010204	20° Tip, Hooked Right	N/A	6.8mm
1010208	20° Tip, Short, Serrated, Left	N/A	6.8mm
1010209	20° Tip, Short, Serrated, Right	N/A	6.8mm



CODE	DESCRIPTION	BITE WIDTH	WORKING DIAMETER
1010210	60° Tip, Hooked Left	N/A	6.8mm
1010212	60° Tip, Hooked Right	N/A	6.8mm



Linear Graspers

CODE	DESCRIPTION	BITE WIDTH	WORKING DIAMETER
1010301	3.4mm, Straight, Cupped	3.4mm	3.7mm



CODE	DESCRIPTION	BITE WIDTH	WORKING DIAMETER
1010303	3.4mm, Straight, Alligator	3.4mm	3.7mm



CODE	DESCRIPTION	BITE WIDTH	WORKING DIAMETER
1010306	3.4mm, Straight, Loose Body	3.4mm	3.7mm

